





LONG-HANDLED TOE WASHER & DRIER

PURPOSE:

To allow independent bathing of feet when reach is restricted.



terry toweling
slip-cover

MATERIALS:

- 1 rubber or plastic bowl scraper or facsimile.
- 1 piece of aluminum (65 ST 6 bar) or equivalent, to lengthen scraper.
- 2 pieces of wood, 3/8" x 1" x 4" for handle.
- 2 pieces of terry toweling 5" x 5" for slip-covers, (1 for washing, 1 for drying), with elastic or velcro around top.

CANADIAN ARTHRITIS SOCIETY